




SKU: 1569007

Rocky 4

The Shopprider Rocky 4 is a robust mobility scooter that features a 4-pole motor, 225kg weight capacity, and front and rear suspension.



Colours



Features



Speed
10 km/h



Weight Capacity
225kg



LED Lighting
Front and rear lighting



Incline Capability
Up to 6° (180kg)
Up to 9° (156 kg)
Up to 12° (136kg)

Full Features

- High back contoured seat with headrest
- Gas adjustable steering tiller
- Loop handle with finger operated controls
- Front and rear LED lighting, including hazard lights
- High / Low speed settings
- Electromagnetic brake
- Front + Rear (Adjustable) suspension
- Large battery capacity
- Front and rear suspension
- Adjustable width and foldable arm-rests
- Sliding / Rotating / Height adjustable seat

Specifications

SPECIFICATIONS	889SL
Overall Dimensions	1430mm (L) x 680mm(W) x 13300mm (H)
Maximum Speed	10km/h
Turning Circle	1300mm
Maximum Range	1200mm
Safe Climbing Angle	9°
Max Climbing Angle	12°
Charger	4 Amp
Motor Size	2.0Hp high torque motor (4 Pole)
Batteries	2 x 75 Amp minimum
Maximum Range	up to 40 km
Wheel Size	330mm - Front + Rear
Turning Radius	1430mm
Product Weight	146 kg (with batteries)
Maximum User Weight	225kg - on flat ground

Warranty

Warranty Information	
Scooter Frames	10 Years Structural Frame 18 Months Drive Train and Batteries
Batteries (Gel/Hybrid)	18 Months
Scooter Parts and Accessories	1 Year

IMPORTANT : Always be sure to read your owner's manual thoroughly before operation.

Silverwater
Unit 8, 110-120 Silverwater Rd
Silverwater NSW 2128
Ph: 9649 2111

Newcastle
7 Poynton Pl
Thornton, NSW 2322
Ph: 4047 3608

activemobility.com.au

