




SHOPRIDER®

Rocky 8

A powerful heavy duty scooter built tough for superb ride quality and smooth handling




Colours Available



Cobalt
Blue



Burgundy
Red



Gotham
Grey

Features



Seat
Large High Back
Contoured Captain
Seat with Headrest



Wheels
4 Wheel



Weight Capacity
225kg
Bariatric Suitable



Lighting
LED Headlights,
Rear Lights, Indicators,
Hazard Lights


Full Features

- Powerful 4 Pole Motor
- Large High Back Contoured Captain's Seat with Headrest
- Gas Adjustable Steering Tiller
- Lift Up & Adjustable Width Arm Rests
- Sliding Rotating Height Adjustable Seat
- Loop Handle with Finger Operated Controls
- Electromagnetic Brake
- Removable Shopping Basket
- Large Battery Capacity
- Horn
- Deluxe LCD Display Dashboard
- Large Spoke Mag Rims
- Front & Rear Adjustable Suspension
- Front & Rear Lights
- Hazard Lights / Indicators
- High / Low Speed Settings
- Exceptional Ground Clearance

Specifications

SPECIFICATIONS	SHOPRIDER ROCKY 8
Length	159cm
Width	73cm
Height	128.5cm
Safe Climbing Angle ¹	12 Degrees
Max Climbing Angle ¹	12 Degrees
Maximum Weight Capacity ¹	225kg
Bariatric Suitable	Yes
Charger	8Amp
Motor Size	2.0Hp (4 pole)
Maximum Speed	10km/h
Maximum Range ²	Up to 54km
Wheel Size	33cm
Turning Radius	186cm
Scooter Weight (with batteries)	147kg

¹Single User



²Driving range varies with the factors shown below:

- (a) The weight of the occupant
- (b) Level ground & ambient temperature
- (c) Battery condition / Type of charger
- (d) Driving Style

IMPORTANT : Always be sure to read your owner's manual thoroughly before operating your machine.

WARNING : Do not exceed 10km/ph on public roads or paths.

Active Mobility
Unit 8, 110-120 Silverwater Rd
Silverwater NSW 2128
PH: 9649 2111
activemobility.com.au



SHOPRIDER