



SHOPRIDER
Viking
 779 XLS

Rugged and durable, the Viking 779 XLS turns heads without compromising on comfort due to a deluxe suspension system



Colours Available



Injection Red

Features



Seat
 High Back
 Contoured Seat
 with Headrest



Weight Capacity
 185kg
 Bariatric Suitable



Wheels
 3 Wheel



Lighting
 LED Headlights with
 low/high beam setting,
 Rear Lights, Indicators,
 Hazard Lights


Full Features

- Powerful 4 Pole Motor
- High Back Contoured Seat with Headrest
- Sliding Rotating Height Adjustable Seat
- 'Slow When Turning' Safety Feature
- Elevated Handlebar with R/H Paddle Control
- LCD Display
- High/ Medium/ Low Speed Settings
- Rear Wheel Drive
- Electromagnetic Brake
- Front Disc Handbrake
- Front & Rear Deluxe Suspension System
- LED Headlight with Low/ High Beam
- Rear Activated Brake Lights
- Indicators
- Hazard Lights / Indicators

Specifications

SPECIFICATIONS	SHOPRIDER VIKING 779 XLS
Length	155cm
Width	70cm
Height	130.5cm
Safe Climbing Angle ¹	9 Degrees
Max Climbing Angle ¹	9 Degrees
Maximum Weight Capacity ¹	185kg
Bariatric Suitable	Yes
Charger	8Amp
Motor Size	2.0Hp High Torque Motor (4 pole)
Maximum Speed	10km/h
Maximum Range ²	Up to 54km
Wheel Size	40.6cm
Turning Radius	165cm
Scooter Weight (with battery)	140kg

¹Single User



²Driving range varies with the factors shown below:

- (a) The weight of the occupant
- (b) Level ground & ambient temperature
- (c) Battery condition / Type of charger
- (d) Driving Style

Warranty Information:

Scooters: 10 year structural frame, 18 months drivetrain and batteries

Parts and Accessories: 12 months

Batteries (Gel/Hybrid): 18 months

IMPORTANT : Always be sure to read your owner's manual thoroughly before operating your machine.

WARNING : Do not exceed 10km/ph on public roads or paths.