



SHOPRIDER
Viking
 779 XLS

Rugged and durable, the Viking 779 XLS turns heads without compromising on comfort due to a deluxe suspension system



Colours Available



Injection Red

Features



Seat
 High Back
 Contoured Seat
 with Headrest



Weight Capacity
 185kg
 Bariatric Suitable



Wheels
 3 Wheel



Lighting
 LED Headlights with
 low/high beam setting,
 Rear Lights, Indicators,
 Hazard Lights


Full Features

- Powerful 4 Pole Motor
- High Back Contoured Seat with Headrest
- Sliding Rotating Height Adjustable Seat
- 'Slow When Turning' Safety Feature
- Elevated Handlebar with R/H Paddle Control
- LCD Display
- High/ Medium/ Low Speed Settings
- Rear Wheel Drive
- Electromagnetic Brake
- Front Disc Handbrake
- Front & Rear Deluxe Suspension System
- LED Headlight with Low/ High Beam
- Rear Activated Brake Lights
- Indicators
- Hazard Lights / Indicators

Specifications

| SPECIFICATIONS | SHOPRIDER VIKING 779 XLS |
|--------------------------------------|----------------------------------|
| Length | 155cm |
| Width | 70cm |
| Height | 130.5cm |
| Safe Climbing Angle ¹ | 9 Degrees |
| Max Climbing Angle ¹ | 9 Degrees |
| Maximum Weight Capacity ¹ | 185kg |
| Bariatric Suitable | Yes |
| Charger | 8Amp |
| Motor Size | 2.0Hp High Torque Motor (4 pole) |
| Maximum Speed | 10km/h |
| Maximum Range ² | Up to 54km |
| Wheel Size | 40.6cm |
| Turning Radius | 165cm |
| Scooter Weight (with battery) | 140kg |

¹Single User



²Driving range varies with the factors shown below:

- (a) The weight of the occupant
- (b) Level ground & ambient temperature
- (c) Battery condition / Type of charger
- (d) Driving Style

IMPORTANT : Always be sure to read your owner's manual thoroughly before operating your machine.

WARNING : Do not exceed 10km/ph on public roads or paths.