



xtaga Gazelle

Power and durability meets exceptional style. Built for adventure, the xtaga Gazelle glides across all terrains with supreme ease and comfort making it the perfect outdoor mobility scooter.



Colours Available



Carmine Red



Champagne Topaz

Features



Seat
High back captain seat with head rest



Weight Capacity
160kg
Bariatric Suitable



Wheels
4 Wheel



Lighting
LED headlight, Indicators, Side & Rear Reflectors, Brake Lights

Full Features

- 4-wheel suspension for a smooth journey over most terrains
- Height and width adjustable comfort padded armrests
- Easy drive adjustable tiller with wrap-around handles
- Contoured, high-back captain seat with padded headrest
- Seat swivel adjustment
- Small turning circle for maneuvering around tight spaces
- Weight capacity up to 160kg
- Front basket
- LED headlight & reflectors
- Reverse warning alarm
- Horn
- Elegant style
- Available in 2 premium colours

Specifications

SPECIFICATIONS	XTAGA GAZELLE
Maximum Weight Capacity	160kg
Maximum Speed ¹	10km/h
Maximum Range ¹	25-35km
Length ²	137cm
Width ²	67cm
Scooter Weight ^{2,3}	75.2kg / 105kg with Battery
Ground Clearance ²	10cm
Seat to Ground	68.5cm (Adjustable)
Turning Radius ²	139cm
Battery Requirements	40Ah
Charger	Powerfirst 5A
Seating	Captain Seat with Headrest
Motor	700W

¹Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition and tyre condition. These specifications can be subject to a variance of (+/- 10%).

² Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or - 3%).

³Battery weight may vary based on manufacturer.

Compliance

COMPLIANCE	
Australian Register of Therapeutic Goods (ARTG)	xtaga Gazelle is entered in the ARTG 406770
ISO	ISO9001, ISO13485

IMPORTANT : Always be sure to read your owner's manual thoroughly before operating your machine.

WARNING : Do not exceed 10km/ph on public roads or paths.

Warranty Information:

Scooters: 3 years (scooter frame)

Parts: 2 years (electronics)

Batteries: 12 Months

Silverwater
Unit 8, 110-120 Silverwater Rd
Silverwater NSW 2128
Ph: 9649 2111

Newcastle
7 Poynton Pl
Thornton, NSW 2322
Ph: 4047 3608

activemobility.com.au

