



xtaga Impala

The ultimate in portability, power and performance. The xtaga Impala is packed with features and perfect for those seeking a compact scooter without having to compromise on luxury or performance.



Colours Available



Puya Teal



Cobalt Blue

Features



Seat
Mid back comfortable seat that can be folded



Weight Capacity
145kg



Wheels
4 Wheel



Lighting
LED headlight, front basket, reflectors, brake lights

Full Features

- Quick disassembly for transporting
- Height and width adjustable comfort padded armrests
- Pneumatic tyres for excellent outdoor performance
- Easy drive adjustable tiller with ergonomic wraparound handles
- Contoured, mid-back captain seat with padded headrest
- Large weight capacity
- Front & rear suspension
- Front basket
- Seat pocket for storage
- LED lower headlight & reflectors
- Horn
- Sleek, sporty style
- Seat swivel adjustment

Specifications

SPECIFICATIONS	XTAGA IMPALA
Maximum Weight Capacity	145kg
Maximum Speed ¹	8km/h
Maximum Range ¹	Up to 29km
Length ²	112cm
Width ²	56cm
Scooter Weight ^{2,3}	59.5kg / 82.5kg with Battery
Ground Clearance ²	9cm
Seat to Ground	47cm (Adjustable)
Turning Radius ²	134.5cm
Battery Requirements	50Ah
Charger	Input AC 100-240V, Output DC24V, 3.5Amp
Seating	Captain Seat with Headrest
Motor	420W

¹Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition and tyre condition. These specifications can be subject to a variance of (+/- 10%).

² Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or - 3%).

³Battery weight may vary based on manufacturer.

Compliance

COMPLIANCE	
Australian Register of Therapeutic Goods (ARTG)	xtaga Impala is entered in the ARTG 406770
ISO	ISO9001, ISO13485

IMPORTANT: Always be sure to read your owner's manual thoroughly before operating your machine.

WARNING: Do not exceed 10km/ph on public roads or paths.

Warranty Information:

Scooters: 3 years (scooter frame)

Parts: 2 years (electronics)

Batteries: 12 Months

Silverwater
Unit 8, 110-120 Silverwater Rd
Silverwater NSW 2128
Ph: 9649 2111

Newcastle
7 Poynton Pl
Thornton, NSW 2322
Ph: 4047 3608

activemobility.com.au

