



# xtaga Impala

The ultimate in portability, power and performance. The xtaga Impala is packed with features and perfect for those seeking a compact scooter without having to compromise on luxury or performance.



## Colours Available



Puya Teal



Cobalt Blue

## Features



**Seat**  
Mid back comfortable seat that can be folded



**Weight Capacity**  
145kg



**Wheels**  
4 Wheel



**Lighting**  
LED headlight, front basket, reflectors, brake lights

# Full Features

- Quick disassembly for transporting
- Height and width adjustable comfort padded armrests
- Pneumatic tyres for excellent outdoor performance
- Easy drive adjustable tiller with ergonomic wraparound handles
- Contoured, mid-back captain seat
- Large weight capacity
- Front & rear suspension
- Front basket
- Seat pocket for storage
- LED lower headlight & reflectors
- Horn
- Sleek, sporty style
- Seat swivel adjustment

# Specifications

SPECIFICATIONS	XTAGA IMPALA
Maximum Weight Capacity	145kg
Maximum Speed <sup>1</sup>	8km/h
Maximum Range <sup>1</sup>	Up to 29km
Length <sup>2</sup>	112cm
Width <sup>2</sup>	56cm
Scooter Weight <sup>2,3</sup>	59.5kg / 82.5kg with Battery
Ground Clearance <sup>2</sup>	9cm
Seat to Ground	47cm (Adjustable)
Turning Radius <sup>2</sup>	134.5cm
Battery Requirements	50Ah
Charger	Input AC 100-240V, Output DC24V, 3.5Amp
Seating	Captain Seat
Motor	420W

<sup>1</sup>Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition and tyre condition. These specifications can be subject to a variance of (+/- 10%).

<sup>2</sup>Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or - 3%).

<sup>3</sup>Battery weight may vary based on manufacturer.

# Compliance

COMPLIANCE	
Australian Register of Therapeutic Goods (ARTG)	xtaga Impala is entered in the ARTG 406770
ISO	ISO9001, ISO13485

**IMPORTANT:** Always be sure to read your owner's manual thoroughly before operating your machine.

**WARNING:** Do not exceed 10km/ph on public roads or paths.

### Warranty Information:

**Scooters: 3 years (scooter frame)**

**Parts: 2 years (electronics)**

**Batteries: 12 Months**

### Silverwater

Unit 8, 110-120 Silverwater Rd  
Silverwater NSW 2128  
Ph: 9649 2111

### Newcastle

7 Poynton Pl  
Thornton, NSW 2322  
Ph: 4047 3608

[activemobility.com.au](http://activemobility.com.au)



**xtaga**