



Azure rollo+

Boasting ergonomic design, supreme adjustability, comfort and class, the rollo+ walker is the perfect combination of function and style in mobility.



Colours Available



Olive Grey

Features



Adjustable Backrest

Up to 80mm of seat depth and 160mm height adjustment



Weight Capacity

180kg



Walker Weight

9.45kg

Full Features

- Suitable for
- Sturdy construction lightweight, aluminum frame weighing under 10kg
- One touch, easy-fold and unfold mechanism for convenient transport and storage
- Safety breaks and push down locks for additional safety and stability when rising from a seated position
- Ergonomically designed hand grips with safety reflectors
- Ultra-comfort seat with supreme cushioning
- Patent adjustable padded backrest providing up to 80mm of seat depth adjustment and 160mm height adjustment
- Height adjustable handles allowing up to 130mm of additional handle height
- Large removable tote bag with zip allows for ample storage of a handbag and belongings
- Easy manouverability and smooth handling over most terrains

Specifications

SPECIFICATIONS	AZURE ROLIO+
Maximum Weight Capacity	180kg
Width	66.5cm
Depth	73cm
Height	92cm
Unit Weight	9.45kg
Folded Width	28cm
Folded Depth	73cm
Folded Height	92cm
Handle Height	92cm - 105cm
Seat Width	56cm
Seat Depth	30.5cm
Seat Height	60cm
Front Wheel	25cm / 10"
Rear Wheel	20cm / 8"

IMPORTANT : Always be sure to read your owner's manual thoroughly before operation

Silverwater
Unit 8, 110-120 Silverwater Rd
Silverwater NSW 2128
Ph: 9649 2111

Newcastle
7 Poynton Pl
Thornton, NSW 2322
Ph: 4047 3608

activemobility.com.au


Active Mobility

rolio+