

Battery Reference Sheet



Name of Chair	Battery Weight* (kg)	Capacity of Battery (Ah)	Type of Battery
Jovy Carbon Wheelchair	25	2 x 24v - 6Ah	Lithium
Shoprider Como	47.4	2 x 12v - 14Ah	SLA Batteries
Streamer	73	2 x 12v - 36Ah	SLA Batteries
Shoprider Venice Light Weight	32	12v - 10Ah	Lithium Iron
Shoprider Puma HD14	120	2 x 12v - 75Ah	SLA Batteries
Puma 10	87	2 x 12v - 36Ah	SLA Batteries
Cougar 10w Power Lift	100	2 x 12v - 50Ah	SLA Batteries
Cougar 10w Power Tilt	100	2 x 12v - 50Ah	SLA Batteries
R40 Fusion	102	2 x 12v - 55Ah	SLA Batteries
Jazzy 623	132	2 x 12v deep-cycle	SLA Batteries
Jazzy Air 2.0	148.56	2 x 12v - 40Ah	SLA Batteries
Jazzy Select 6	108	2 x 12v - 30-34Ah	SLA Batteries
Jazzy Select Elite HD	121	2 x 12v - NF - 22Ah	SLA Batteries
Pride Go Chair	71.72	2 x 12v - 12Ah	SLA Batteries
Shoprider Aspire Mini	100	2 x 12v - 26Ah	SLA Batteries
Shoprider Aspire Maxi	126	2 x 12v - 40Ah	SLA Batteries
Vector	147	2 x 12v - 55Ah	SLA Batteries
M12 MWD Power Wheelchair	106	2 x 12v - 55-62Ah	SLA Batteries
Glide Series 4	52	2 x 12v - 24Ah	SLA Batteries
Glide Centro	110	2 x 12v - 75Ah	SLA Batteries
Quickie Q400 Plus	130	2 x 12v - 60-72Ah	SLA Batteries
Centro- Vertical Lift	110	2 x 12v - 70-90Ah	SLA Batteries
Centro- Lift- Anterior Reach	158-168	2 x 12v - 70-90Ah	SLA Batteries
Glide Centro Bariatric Power	134-142	2 x 12v - 75Ah	SLA Batteries
Glide Centro XT	110	2 x 12v - 76Ah	SLA Batteries
Glide Series 7	110	2 x 12v - 36Ah	SLA Batteries
All Track HD3	90	2 x 12v - 70Ah	SLA Batteries
All Track M3	94	2 x 12v - 70Ah	SLA Batteries
All Track R3 Hybrid	91	2 x 12v - 50-70Ah	SLA Batteries
Sango Advanced	120	2 x 12v - 60Ah	SLA Batteries
Sango XXL	120	12v - 78-100Ah	SLA Batteries
Sango Slimline	120	2 x 12v - 60Ah	SLA Batteries
Q500 M Sedeo Pro	130	2 x 12v - 61Ah	SLA Batteries
Zippie ZM-310 Kids	80	2 x 12v - 55Ah	SLA Batteries
Q700 M Sedeo Ergo	140	12v - 60/72Ah	SLA Batteries
R-44	100	2 x 12v - 70Ah	SLA Batteries
Jazzy 1450	160	2 x 12v - 70Ah	SLA Batteries
Q610	80	2 x 12v - deep cycle	SLA Batteries
Q6 Edge HD	100	2 x 12v 70Ah	SLA Batteries
Quantum 1450	130	2 x 12v - 70Ah	SLA Batteries
Q6 Edge 2.0	105	2 x 12v - 70Ah	SLA Batteries
TA IQ MWD	110	2 x 12v - 70Ah	SLA Batteries
TA IQ MWD	110	2 x 12v - 70Ah	SLA Batteries
TA IQ MWD	110	2 x 12v - 70Ah	SLA Batteries

* Weights are a reference only. Accessories and functions can be fitted to the chairs and therefore increase the weight